Develops individual finger strength, flexibility and coordination
Also builds hand and forearm strength

If the individual is being treated for a hand disorder, specific exercises and the intensity and number of repetitions should be approved by a healthcare professional.

Following injury or surgery, consult a healthcare professional before performing exercises.

If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.

Exercises should be done at a slow, comfortable pace. Suggested exercises follow.

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Art.Nr. 2556606 ff.


FLEX-ION

Approved by Certified Hand Therapists

FLEX-ION

Approved by Certified Hand Therapists

Resistance level: per finger / per hand	extra light	0,35 - 1,1 kg	0.75 - 2.5 lbs
Widerstandskraft: pro finger / gesamtes	light	0,7 - 2,3 kg	1.5 - 5.0 lbs
Niveau de résistance: par doigt / sur l'ensemble	medium	1,4 - 4,5 kg	3.0 - 10.0 lbs
Weerstandsniveau: per vinger / geheel	heavy	2,3 - 7,3 kg	5.0 - 16.0 lbs
Nivel de resistencia: por dedo / sistema completa	extra heavy	3,2 - 10,4 kg	7.0 - 23.0 lbs
Livello di resistenza: per singolo gito / intero apparecchio	super heavy	4,1 - 14,1 kg	9.0 - 31.0 lbs

Daily use conditions and rehabs the hands, wrists and forearms by flexing and toning the muscles, ligaments and tendons. Flex-ion provides comprehensive rehabilitation to regain healthy hands.

EXERCISE EACH FINGER INDIVIDUALLY, OR THE ENTIRE HAND

For additional exercises, developed by a licensed hand therapist, please visit www.msdband.com

MSD Europe
Westdijk 150
2830 Tiselt
Belgium
www.msdband.com
info@msdband.com
Find us on

Contains one unit
Exercise program enclosed

MSD North America, Inc.
1100 Grape Avenue
St. Cloud, FL 34769
United States of America
www.msdband.com
info@msdband.com

**EN - Hook Position**

Place finger tips on individual buttons and position ergonomic palm bar with the hook over the web space between the thumb and index finger. Keeping finger tips slightly bent, flex all fingers toward the center of the unit.

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EN - Tip-to- Tip Pinch

Place the tip of the thumb on the center of the ergonomic palm bar. Position the finger to be exercised on one of the center buttons. Keeping both tips slightly bent, flex the thumb and the finger being exercised toward the center of the unit at the same time.

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**EN - Gross Grasp**

Place the middle of the fingers on the buttons with the finger tips wrapped over the top. Positioning the hook of the ergonomic palm bar over the web space between the thumb and index finger, flex the thumb and all fingers toward the center of the unit as if making a fist.

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EN - Fingertip Flexion

Place all finger tips on their own button, with the tip of the thumb on the center of the ergonomic base. Keeping all finger tips bent, flex the thumb and all the fingers toward the center of the unit at the same time.

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**EN - Key Pinch**

Place the hook of the ergonomic palm bar on the inside of the bent index finger. With the tip of the thumb slightly bent, flex the thumb on the first button.

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For additional exercises, developed by a licensed hand therapist, please visit www.msd-band.com

EN - Trigger Pinch

Loosely hold the unit vertically with the buttons toward the fingers and the hook of the ergonomic palm bar positioned over the web space between the thumb and index finger. With the middle of the index finger, press the second button as if squeezing a trigger.

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