SISSEL® INTENSIVE MASSAGER

GEBRAUCHSANWEISUNG - bitte aufmerksam lesen und gut aufbewahren
OPERATING MANUAL - please read carefully and save
MANUEL D'UTILISATION - lisez attentivement ce manuel et conservez le
MANUAL DE INSTRUCCIONES - lea con atención las instrucciones y guárdelas
ISTRUZIONI PER L'USO - leggere attentamente e conservare con cura
GEBRUIKSAANWIJZING - aandachtig lezen en bewaren
BRUKSANVISNING - läs denna noga och förvara den på ett säkert ställe
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - Παρακαλούμε διαβάστε προσεκτικά και κρατήστε το παρόν
έγγραφο

SISSEL® Intensive Massager

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SISSEL[®] Intensive Massager

Congratulations on your purchase of the latest development in hand held massagers. Developed by health professionals, the SISSEL[®] *Intensive Massager* was designed to assist in the therapeutic relief of muscle tension at home.

The device is not intended for professional use. It is not a replacement for proper medical care. Please consult your physician for any pain of unknown cause.

Thank you for choosing SISSEL[®] *Intensive Massager*.

Important Safety Instructions

When using an electrical appliance like SISSEL[®] *Intensive Massager*, basic precautions should always be followed, including the following:

1. Always unplug the SISSEL[®] *Intensive Massager* from the electrical outlet immediately after using and before cleaning
2. Do not reach for an appliance that has fallen into water. Unplug immediately.
3. Do not use while bathing or in a shower.
4. Do not place or store the SISSEL[®] *Intensive Massager* where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.

Warning:

To reduce the risk of burns, fire, electric shock or injury to persons:

1. SISSEL[®] *Intensive Massager* should never be left unattended when plugged in. Unplug SISSEL[®] *Intensive Massager* from outlet when not in use, before putting on or taking off parts, and before cleaning.
2. Do not operate under a blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injuries to persons.
3. Close supervision by adults is necessary when SISSEL[®] *Intensive Massager* is used by, on, or near children. Treatment on invalids, and mentally and sensory impaired or disabled persons only after consulting a physician. Do not use on young children (0-10 years).
4. Use SISSEL[®] *Intensive Massager* only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
5. Never operate SISSEL[®] *Intensive Massager* if it has a damaged cord or plug, if it is not working properly, if it has been dropped, damaged, or dropped into water. Return the appliance to a service center for examination and repair.
6. Do not carry SISSEL[®] *Intensive Massager* by supply cord or use cord as a handle.
7. Keep the cord away from heated surfaces.
8. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
9. Never operate when ventilation openings are blocked or covered with lint, hair etc.

Power supply: SISSEL[®] *Intensive Massager* operates with ordinary 230 volt household current. It is too powerful to run on batteries and is double-insulated, it has two systems of insulation instead of grounding. Therefore maintenance and repair shall only be done by qualified personnel.

Care / Maintenance: If the external surfaces become soiled, wipe clean with a cloth dampened with a mild soap. Do not use thinner, benzene, or any other solvents. SISSEL[®] *Intensive Massager* is maintenance free, constructed with permanently lubricated sealed ball bearings. No additional lubrication is necessary. SISSEL[®] *Intensive Massager* has no user serviceable parts inside. **Opening by anyone other than an authorized service center will void the warranty.**

We recommend using the SISSEL[®] *Intensive Massager* for no longer than 20 minutes at one time, then allowing at least 20 minutes for cooling. Overuse, as well as applying pressure, could cause the product to overheat. Should this occur, discontinue use and allow the unit to cool before operating.

SISSEL[®] *Intensive Massager* has an automatic reset thermal switch. If this appliance overheats, the thermal switch will trip and it will cease operating. Shut power off, allow sufficient time for cooling and the appliance with restore itself to normal operation.

When not in use, store the cord coiled with appliance. When in use, ensure the cord is untangled, free and not twisted. For storage, loop the cord several times into 30cm loops and tightly wrap the cord tie around the bundle.

Important: Massage should be pleasant and comfortable and especially so with SISSEL[®] *Intensive Massager*. Its percussive action produces a feeling of steady, thorough massage reaching down through your muscles. If you have any concerns regarding your health, consult your doctor before using SISSEL[®] *Intensive Massager*.

Individuals with pacemakers or in case of pregnancy, diabetes or illness, consult your doctor before using SISSEL[®] *Intensive Massager*.

Itching may occur if:

The device is left in one spot too long (over 60 seconds, due to increased circulation bringing warm blood to the area, if the user has poor circulation or if irritating material or clothing (e.g. synthetic materials) is used between the Intensive Massager and the skin.

SISSEL[®] *Intensive Massager* should not be used over:

- Abdomen or low back during pregnancy.
- Swollen or inflamed areas or skin eruptions (bruises, cuts, fractures, etc.)
- Testicles, ovaries, breast tissue, eyes, ears, head or neck
- Directly onto bone
- Varicose veins
- Do not use on unexplained calf pain.

If you feel pain or discomfort, discontinue use immediately and consult your physician.

HOW TO USE

- SISSEL[®] *Intensive Massager*'s percussive action, delivers a thorough and effective deep massage. Its pulsating energy penetrates even large body parts with ease. When being applied straight down (i.e. using gravity) its own 1.4 kg. weight supplies all the required force. Adding extra weight is not recommended and may cause overheating.
- SISSEL[®] *Intensive Massager* is designed for massage of soft tissue (muscle) not hard tissue (bone) and goes after muscle tension with a percussive stroke that travels downward into your muscles, releasing muscular tension. The muscle SISSEL[®] *Intensive Massager* massages should be at least 1/8" thick or SISSEL[®] *Intensive Massager* may feel too strong and bounce or shake. Larger, thicker muscles are best with SISSEL[®] *Intensive Massager*. For thinner areas, turn down the power setting to a comfortable level

- If stroking the legs, always stroke toward the heart. Do not use on unexplained calf pain or varicose veins if you experience pain in a muscle or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.

Attachments: SISSEL® *Intensive Massager* has 2 sets of interchangeable massage heads, with varying densities, to allow you to massage all areas effectively. Select the 'hardness' and 'shape' of the massage attachment that you prefer. These heads snap on and off

Power Settings: SISSEL® *Intensive Massager* has a rotary dial for variation of the power/speed settings. Select the speed you prefer, for the muscle you wish to release. If it feels 'too strong', lower the speed setting to a more acceptable level.

Strength settings:

- High Strength: Dial power setting #5 and/or using the black massage heads.
- Medium Strength: Dial power setting #3.
- Low Strength: Dial power setting #1 and/or using the gray, massage heads.
- Extra Low Strength: Massaging through a folded towel will reduce its impact even below the lowest power setting.

Product features:

- Variable Speed Control. Dial allows adjustment of power from 1,000 to 3,000 pulses per minute.
- Extended handle to reach down your back.
- 2 sets of Interchangeable Massage Heads. Added flexibility of choosing the firmness of massage that is just right for you, for your individual massage preference. The heads simply snap on and off with a slight twisting motion.

1. Soft rounded for gentle massage.

2. Firm surface for harder massage

- NOTE: DO NOT USE SISSEL® *Intensive Massager* WITHOUT A SET OF MASSAGE HEADS.
- Easy to use. Just rest it on your tight muscles.

Product specifications:

- Weight: 1.4kg
- Volts: ~ 220-230 Volt/ 50 Hz
- Cable length: 3.7m
- Pulses: 1,000 to 3,000 per minute
- Motor: DC with sealed ball bearings
- Protection class: II

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SISSEL® GmbH