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Italia

Country

HEREWITH DECLARES, UNDER SOLE RESPONSIBILITY, THAT FOLLOWING MEDICAL DEVICE:

Item name	HORIZONTAL LEG PRESS
Item code	MCC.07/KGG
Line	Ext3Pro / KGG Set
Class	Counterweighted equipment
Planned use	Horizontal Leg Press for physiotherapeutic exercises

**IS COMPLIANT TO ALL RELEVANT INSTRUCTIONS OF REGULATION 2017/745/EU
REGARDING MEDICAL DEVICES**

It has been applied the attached VIII of Regulation 2017/745/EU to determine the classification of medical device, which results in **Class I**.

REGULATIONS AND COMPLIANCE

- UNI EN 913:2021 Gymnastic equipment - General safety requirements and test methods.
- UNI EN ISO 20957-1:2013 General safety requirements and test methods for stationary training equipment.
- UNI EN ISO 20957-2:2021 Stationary training equipments – strenght training equipments, specific safety additional requirements and testing methods.
- UNI EN ISO 20957-4:2017 Strenght training benches, additional specific safety requirements and test methods.
- We also declares that:
 - The year of production is reported on the label (with CE mark) put on the equipment.
 - The authorized person to guard the technical file of construction at our plant, following what is expected by the Regulation, is Terziotti Rubes as legal representative.
 - The serial number of the equipments is reported on the same label.

Pegognaga – Italy, May 30th 2024
Place and date

TORMEC SRL
Via G. Rossa, 9/11 - 46020 Pegognaga (MN)
Tel. 0376/558600 - Fax 0376/559995
C.F. P.IVA e Reg. Impr. MN 01661260206
R.E.A. n. 173331
Company stamp and legal representative's signature

Sport-Thieme

Medical functional trainer

Thank you for choosing this product!

Please read the following notes carefully before use so you can enjoy this product safely and keep these instructions for future reference. If you have any questions, our team is here for you.

Introduction

Thank you to have chosen Sport-Thieme Medical functional trainer. All the strength machines are projected and built to satisfy all the users' needs for a handy, correct and safe use.

This user manual contains little suggestion to use the device and precious advices for your safety.

It is suggested to read carefully the entire manual before using the medical device. In case of any doubts, please contact the distributor, who is able to help you.

Avilability of the manual:

User manual can be required at sport-thieme.de/service/bedienungsanleitungen it can be requested when the product is over its useful life. The instructions are attached to the devices and they accompany them for their entire useful life.

Assembly

About Product

Technical features:

Equipment weight :	350 kg
User's max weight:	130 kg
Max load:	150 kg (0-20 kg: plates 5/5 kg ; 20-150: plates 10/10 kg).
Assembled equipment dimensions :	length 220 cm x width 120 cm x height 200 cm
Material:	Steel frame

General description

Sport-Thieme Medical functional trainer allows a comfortable and practical movement of legs stretching, involving quadriceps, glutei, rectus femoris and hip adductors in most accentuated flexions. The seat with release device permits a comfortable starting position with extended legs. Handles allow a steady positioning. The big footrest permits permits to manage the inclination of feet resting. The flowing of carriage on bearings makes the effort extremely fluid and equilibrated.

Field of application / Installation condition

The use of Sport-Thieme Medical functional trainer to indoor use in professional healthcare facilities or in the patient's home.

- Clinics (patient rooms, intensive care),
- Physiotherapy

Field of application / Installation condition

The use of Sport-Thieme Medical functional trainer to indoor use in professional healthcare facilities or in the patient's home.

- Clinics (patient rooms, intensive care),
- Physiotherapy clinics. Furthermore:
- Do not use Sport-Thieme Medical functional trainer in wet, humid or hot environments.
- Sport-Thieme Medical functional trainer only on a flat, steady and not slippery surface.
- Do not use Sport-Thieme Medical functional trainer in open spaces or outdoor areas.

Unpacking

Unpack Sport-Thieme Medical functional trainer:

- Take Sport-Thieme Medical functional trainer out of the packing.
- Check Sport-Thieme Medical functional trainer and all the components /accessories in order to verify that device has not been damaged during transport.
- Check that supply is complete.
- Please immediately communicate any occurred damage to supplier and carrier.

Assembly

Before installing or using the equipment, please wear personal protective equipment, as follows:

- Non slip-shoes
- Protection gloves
- Anti-entanglement suit


Make sure that appointed workers for installation and assembly are equipped with all the personal protective equipment indicated.

Make these actions with help of at least 4 operators, paying attention to avoid physical injuries and/or damages to the equipment while moving the parts.

The device Sport-Thieme Medical functional trainer is supplied in two parts (BASE with sliding seat and WEIGHT STACK), packed with a plastic heat shrinkable cover and fixed to pallet (2x).

Open the packing and lower down the base and the weight stack, on a perfectly flat surface in a suitable position to fix each other through the connection pipes.



Place frontal and posterior connecting pipes according to corresponding holes (marked with a number or a letter) and proceed with fixing by specific screws and washers (already placed into holes of connecting pipes).

Assembly

After having fixed base and weight stack, it needs:

- unroll traction belt and placing it following this sequence:
 - firstly on inferior pulley,
 - then on upper pulley (above weight stack),
 - and lower it on chromed fixing support on weight stack.
- On the weight stack, proceed to fix the end of the belt using the specific spring pin, looking at the following pictures (in sequence):

After the end of assembly, take of the white film that covers front plexiglass protection screens.

Spare Parts And Accessories

Available spare parts for Sport-Thieme Medical functional trainer are the following:

Reference	Quantity	Description
1	2x	guide bushing
2	2x	rod for sliding
3	1x	upper pulley
4	12x	load plate – kg. 10
5	1x	knob for weight selection
6	2x	guides for carriage sliding
7	1x	lower wheel (footrest side)
8	5,3x	belt (linear meters)
9	2x front, 1x rück	protection plexiglass screen
10	1x	lower wheel (weight stack side)
11	1x	intermediate relaunch pulley

Maintenance

Cleaning and disinfection

Cleaning: use only a humid cloth and neutral soap, then dry well before use. Use water warmed by a temperature not over 30°C.

Disinfection: if necessary make a disinfection of the product using a common disinfectant cleaner.

N.B. Absolutely avoid using acidic, alkaline products or solvents such as acetone or thinner.

For patient's and supervisor's safety it is recommended (at least once at year) a suitability check of the device.

Periodic inspection means an examination of the following parts:

- Frame of the device ;
- Check all the welded points ;
- Check fixing of all screws and nuts ;
- Check of integrity of seat, handles and knobs. Furthermore, check gas springs, plexiglass protections, belts and pulleys.

In case of repair, only original parts and accessories must be used.

More informations

Warranty

Excluded from our warranty are commercial buyers who purchase the goods for the purpose of resale.

We generally grant a 3-year warranty, unless stated otherwise on the item. Your statutory warranty rights, the exercise of which is free of charge, remain unaffected by our warranty. Our warranty applies for the period from the handover of the goods to the buyer until the expiration of the respective warranty period indicated. In the event of a warranty claim, the buyer may send the goods to the following address at our expense. At our discretion, the buyer will receive a repair or a replacement delivery. The warranty does not cover normal wear and tear or damage caused by the buyer.

The guarantor is Sport-Thieme GmbH, Helmstedter Straße 40, 38368 Grasleben, Germany.

Warning signs/Notices

Users

Patient

- People who does a first training or recurring training using Sport-Thieme Medical functional trainer.
- Patients are trained by specialised personnel.
- Patients can not train other people in using Sport-Thieme Medical functional trainer.

Specialised Personnel

- People who are skilled through a specialised training in medical or technical-medical sector (for example doctors, therapists)
- People who has professional experience and they did specialised training, therefore they are able to evaluate safety and dangers about use of this kind of medical device.
- People trained by the producer (for example distributors) who are able to supply instructions about the use of the device.

Monitoring

Workouts are made under monitoring by a trained specialist.

Workout under monitoring:

- Specialised personnel monitors the patient during entire workout.
- In this way, specialised personnel can intervene immediately, in case some dangerous situations occur.

Workout without monitoring:

- Specialised personnel stays in the room during entire workout (not directly monitoring the patient).
- Trained specialised personnel can evaluate if the entire workout is made in a safe way.
- Specialised personnel can intervene immediately, in case some dangerous situations occur.
- The use of Sport-Thieme Medical functional trainer is allowed, without monitoring, only if:
 - Patient is able to stand up and keep the position in active way.
 - It is certain that during workout there are no risk of losing consciousness.
 - Patient is able to activate Sport-Thieme Medical functional trainer.

End of Use

Medical devices Sport-Thieme Medical functional trainer are intended to be used in short duration physiotherapy treatments ; they promote the improvement of health through specific physical exercises physiotherapeutic for both muscles and joints.

- A different kind of use of the device (not described in this manual) is forbidden.
- Tormec srl declines any responsibilities for damages caused by improper use of the device or a use other than the one indicated on this manual.
- The producer reserves the right to make changes to the device and the user manual without prior advice, on purpose of improving its features.
- In order to solve pathologies it is necessary that strength machines are used under monitoring of a physiatrist or a corresponding specialist.
- Strength machines are projected and built to be used for promoting improvement of health through specific physical exercises physiotherapeutic for both muscles and joints.

Warning signs/Notices

Medical use

Sport-Thieme Medical functional trainer is a non-motorized training device, which provides a safe environment for training the lower limbs and balance (postural control) and generally retraining muscles and joints.

Indications and target groups of patients

Sport-Thieme Medical functional trainer is suitable to be used mainly in the following situations:
Sport-Thieme Medical functional trainer it is also indicated to allow a muscular and joints rehabilitation, for example to prepare the training of deambulation or to obtain positive effects in improving movements of body.

About therapy by Sport-Thieme Medical functional trainer are mainly applied the same indications and contraindications valid for manual therapeutical treatments. In order to avoid any risks for patient, it is essential to know contraindications. Before any therapeutical session, please check about contraindications that can be occurred.

Contraindications

Sport-Thieme Medical functional trainer must not used by:

- Patients with muscle contractures (deficit of hip flexion / extension and knee > 30°)
- Patients with open wounds that would come in contact with the device (for example, ulcers).
- Lack of stability of bones / pathologies of bone metabolism (for example: fractures, osteoporosis, osteogenesis imperfecta).
- Patients with acute inflammation on lower limbs.
- Patients with acute thrombosis.

Please use Sport-Thieme Medical functional trainer after consulting a physician in the following cases:

- Acute pain during movement.
- Epilepsy.
- Hip dislocation.
- Cardiovascular diseases.
- Patients with acute neurological disease.

Side effects

In infrequent cases:

- Pain (Load pains).
- Skin damages (for exaple: splitting of the skin against the device).

Not compliant use

Sport-Thieme Medical functional trainer is not compliant to do:

- Medical diagnosis.
- Monitorings.
- Measurements.

Do not use Sport-Thieme Medical functional trainer

- Into environments with not compliant conditions to requirements.
- In not monitored workouts, if patient is not able to activate the device by itself.

Biocompatibility

All the components and accessories of Sport-Thieme Medical functional trainer, which patient is constantly in contact with during the expected use, are projected in a way to be compliant to requirements of biocompatibility provider for by currents law.

Warning signs/Notices

Obligation to report unwanted incidents

All serious incidents which happen in relation to the device must be immediately reported to the producer and to competent authority.

Remaining Risks And Safety Stickers

According to current state of knowledges, there are not known and foreseeable risks, neither eventual undesired side effects. Therefore remaining risks are acceptable than evaluated benefits for patient / user coming from performance of the device in normal conditions of use. In order to keep under control risks connected to use of medical devices, workouts must be done with surveillance by a trained specialist.

Workout under surveillance:

During the entire workout, specialised trained personnel stays beside the patient.

- In this way the specialised trained personnel can immediately intervene, in case during workout some dangerous situations occur.

Workout without surveillance:

- During the entire workout, specialised personnel stay in the room (not directly close to the patient).
- Specialised personnel can evaluate if the entire workout is executed in a safe way.

It is believed that remaining risks of medical device are few, most of all considering their classification (all the devices are in Class 1). This because no one of this kind of medical devices can be considered: invasive, active or implantable.

On strength machines equipped with weight stack there is crushing hazard due to plates in movement. To avoid this kind of danger, weight stacks have plexiglass screens for protection.

List of pictograms on the device:

Please check that all the equipments are provided of this pictogram about respecting contents of user manual and instructions.

Please check that strength machines are provided of this pictogram about danger of crushing due to plates in movement.

Warning signs/Notices

Safety

General informations

- Read and observe the user manual.
- Missed reading or lack of compliance to the user manual and maintenance instructions can cause a potential dangerous situation for patients and health workers, with consequent injuries and damages to location.
- Before using, please ask informations to an authorized supplier.
- The assistant must consult specialised personell, taught by the producer.
- Keep an appropriate distance from walls and other obstacles.
- For a safety and effective use, please use Sport-Thieme Medical functional trainer only in the way indicated by following instructions and installation conditions required.
- Please verify thar Sport-Thieme Medical functional trainer is protected from any fluids penetration.
- Do not leave the medical device to unattended children.

Before workout

Before beginning the workout, please check that Sport-Thieme Medical functional trainer works properly. Furthermore:

- Before beginning any workout, please consult a physician to clear if training by Sport-Thieme Medical functional trainer is suitable for patients with cardiovascular disease.
- Before beginning a workout, please agree with a physician the suitable workout parameters for the patient (avoid excessive stress).
 - Adjust mobility according to patient's therapeutic goals.
 - Patient must not do excessive efforts.
 - If patient do not achieve the range of mobility, please change the settings of the device.
- Make sure that no body parts get caught in the equipment while adjusting the Sport-Thieme Medical functional trainer

After workout

- Clean the device using only a humid towel and neutral soap, then dry well before use.
- If necessary, make a disinfection of the product using a common disinfectant cleanser.

Use of device

It is possible to adjust inclination of backrest, acting on the handle positioned on the right of the seat.

It is possible to adjust starting position (and then the extension of legs) acting on handle positioned on the left of the seat.


It is possible to adjust load through the orange knob with selection pin, putting it into the hole on the plate marked with desired weight.

It is possible to adjust inclination of footrest by the specific orange knob on frontal adjusting rod (look at the picture).

Use of device

Muscles involved:

- 1) Adductors
- 2) Gastrocnemius
- 3) Gluteus maximus
- 4) Fibularis longus
- 5) Semimembranosus
- 6) Semitendinosus
- 7) Triceps surae
- 8) Tensor fasciae latae

Moviments and positions:

1. Adjust the backrest, keeping a right positions for hips.
2. Keep the feet in line with hips and keep hands on handles. Keep head, shoulders and back on backrest.
3. Push the low part of the back on backrest, bend the kness and go down until thigh and tibia make a 90° angle.
4. Do not pull up glutei from the seat during the movement.

Safety rules during use:

- Use a proportional load, do not exceed.
- Avoid excessive movements of hip rotation.
- Do not extend completely the knees.
- Keep heels on footrest.
- Steady the return, avoid sudden movements.
- Avoid any kind of contact (user and other people as well) with weight stack during exercise.

TORMEC SRL
Via Guido Rossa 9/11
46020 Pegognaga (MN) - Italy
Tel. 0376/558600 Fax 0376/559995
www.tormec.info e-mail: tormec@tormec.info
P.IVA , C.F. e Reg.Impr.MN n° 01661260206

DECLARATION OF CONFORMITY

(Attached IV Regulation 2017/745/EU)

THE MANUFACTURER

Tormec s.r.l.

Company

Via Guido Rossa, 9/11

Address

Mantova

Province

46020

ZIP code

Pegognaga

City

Italia

Country

HEREWITH DECLARES, UNDER SOLE RESPONSIBILITY, THAT FOLLOWING MEDICAL DEVICE:

Item name	HORIZONTAL LEG PRESS
Item code	MCC.07/KGG
Line	Ext3Pro / KGG Set
Class	Counterweighted equipment
Planned use	Horizontal Leg Press for physiotherapeutic exercises

IS COMPLIANT TO ALL RELEVANT INSTRUCTIONS OF REGULATION 2017/745/EU
REGARDING MEDICAL DEVICES

It has been applied the attached VIII of Regulation 2017/745/EU to determine the classification of medical device, which results in **Class I**.

REGULATIONS AND COMPLIANCE

- UNI EN 913:2021 Gymnastic equipment - General safety requirements and test methods.
- UNI EN ISO 20957-1:2013 General safety requirements and test methods for stationary training equipment.
- UNI EN ISO 20957-2:2021 Stationary training equipments – strenght training equipments, specific safety additional requirements and testing methods.
- UNI EN ISO 20957-4:2017 Strenght training benches, additional specific safety requirements and test methods.
- We also declares that:
 - The year of production is reported on the label (with CE mark) put on the equipment.
 - The authorized person to guard the technical file of construction at our plant, following what is expected by the Regulation, is Terziotti Rubes as legal representative.
 - The serial number of the equipments is reported on the same label.

Pegognaga – Italy, May 30th 2024
Place and date

TORMEC SRL
Via G. Rossa, 9/11 - 46020 Pegognaga (MN)
Tel. 0376/558600 - Fax 0376/559995
C.F. P.IVA e Reg. Impr. MN 01661260206
R.E.A. n. 173331
Company stamp and legal representative's signature

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REGULATIONS AND COMPLIANCE

- UNI EN 913:2021 Gymnastic equipment - General safety requirements and test methods.
- UNI EN ISO 20957-1:2013 General safety requirements and test methods for stationary training equipment.
- UNI EN ISO 20957-2:2021 Stationary training equipments – strenght training equipments, specific safety additional requirements and testing methods.
- UNI EN ISO 20957-4:2017 Strenght training benches, additional specific safety requirements and test methods.

• We also declares that:

- The year of production is reported on the label (with CE mark) put on the equipment.
- The authorized person to guard the technical file of construction at our plant, following what is expected by the Regulation, is Terziotti Rubes as legal representative.
- The serial number of the equipments is reported on the same label.

Pegognaga – Italy, May 30th 2024

Place and date

TORMEC SRL
Via G. Rossa, 9/11 - 46020 Pegognaga (MN)
Tel. 0376/558600 - Fax 0376/559995
C.F. P.IVA e Reg. Impr. MN 01661260206
R.E.A. n. 173331

Company stamp and legal representative's signature